


Model 2435 PSA

AMEL
INSTRUMENTS



The pollution from heavy metals is becoming a crucial point between the environmental dangers derived by production processes. An easy to run and sensitive analysis method is then mandatory for both surveillance and routine controls by the plant managers and pollution control authorities.

Our Model 2435 Tracescan is the solution to all these tasks, with a fully automatic measuring cycle and a straight connection between the instrument and the host PC.

Between the key performances of the instrument you find:

- detection capability down to ppb levels
- no sample digestion or preparation required in most cases
- extremely fast results
- compact instrument with optional battery for in-field operations
- The software given with the instrument allow simple and fast operation from both skilled and unskilled personnel.

Check our wide detection capability library or ask for your specific problem.

From instrument to a PC in few seconds!



Model 2435

Technical specifications

Field of application

PSA technique can be successfully employed in both industrial and research field
Specific analysis of pollutant metals such as: Pb, Cd, Zn, Cu, As can be carried out on

Food and beverages	wine, juice and milk
Water Treatment	tap water, fresh water, sea water and waste water
Galvanic industries	galvanic bath
Pharmaceutical industries	raw material using for packaging
Research	Windows software supplied with Tracescan is completely open: researcher can built his own analysis simply filling-up set up parameter menu

Control Unit

Three-electrodes system

Output voltage compliance	±12V
Response time	≥100µs
Output impedance	≥100Mohm
Ranges	±100µA, ±10µA, ±1µA, ±100nA
Resolution	1/4096 of the selected full scale value

Electrometer

Input voltage	±10V
Input impedance	≥10 ¹² ohm
Input leakage current	≤10pA
Equivalent input noise	≤25µVpp

Analog-to-digital converter

Fast A/D conversion

Potential range	±4.096mV
Resolution	±4.096mV
Conversion rate	≥50ksample/s

Analytical cell

Dosing system	two-way
Stirrer	computer controlled propeller type stirrer
Electrodes	6mm glass body standard N6 conical joint
Reference	Ag/AgCl
Counter-electrode	platinum wire
Working	Glassy carbon 3mm disk (gold and platinum also available)

Techniques implemented

PSA	Potentiometric stripping Analysis
CCSA	Constant current Stripping Analysis

General specifications

Power supply	105, 240V.A.C. 50/60, Hz 30VA (12 Vcc with battery option fit)
Dimensions	300x300x360 mm
Weight	3 kg (5 kg with battery option)

Operating principles

The Potentiometric Stripping Analysis joins the advantages of electrochemistry techniques, like low cost and fast response, to outstanding performance in the interference rejection from unwanted elements. This allows trace analysis directly into the practical sample, without the need for mineralization or sample preparation (as long as the sample is in liquid form).

The technique is based on an accumulation step, identical to the classical voltammetric technique, and a stripping phase that is performed exclusively by a chemical reaction, that is the stripping off the absorbed species in the working electrode by the action of a strong oxidizing solution.

By monitoring the working electrode, the different metals are read in sequence at their own electrochemical potentials: this process is quite insensitive to almost any interfering organic matter and dissolved oxygen too, overcoming traditional voltammetric analysis limitations.

The concentration of the given component is then proportional to the time by current product at the specific redox potential.

The software suite takes care completely of the experiment's execution and data analysis, while the automatic syringe avoids the need of any other equipment for obtaining clean and fast results.