

Lead

Function: Differential Pulse Stripping Voltammetry (DPS/a)

Start Potential (mV)	-800
End Potential (mV)	-200
Current range	1,024 μ A
Scan Speed (mV/s)	30
Deposition time (s)	120
Deposition Pot. (mV)	-800
Number of cycles	3
Delay before sweep (s)	5
Purge and stir time (s)	20
Stirring speed (rpm)	300
Drop Size (a.u.)	60

Lead concentrated standard solution (1 g/l)

Dissolve 1.5986 g of $\text{Pb}(\text{NO}_3)_2$ (pure and dried), in 1 l of 1 % HCl, in a volumetric flask. ($\text{MM}_{\text{Pb}(\text{NO}_3)_2} = 331.21$; $\text{MM}_{\text{Pb}} = 207.2$).

Supporting Electrolyte

A- 37% HCl

B- 1 M $\text{H}_2\text{C}_2\text{O}_4$ and 2 M HCl solution

Dissolve 90 g of $\text{H}_2\text{C}_2\text{O}_4$ (or 126 g of $\text{H}_2\text{C}_2\text{O}_4 \cdot \text{H}_2\text{O}$) and 167 ml of 37% HCl in 1 l of distilled water. Store in a polythene bottle.

Procedure

Add 20 μ l of 37 % HCl to 10 ml of neutralised sample.

Alternatively: add 1 ml of B solution (especially if copper has to be analysed in the same solution).

Analyse sea water, high salt content sample and acidic solution (at pH between 1 and 3) avoiding the addition of the supporting electrolyte.

Samples at pH above 7 are to be neutralised before the addition of the supporting electrolyte.

Working standard solution (1 mg/l)

Dilute 100 μ l of 1 g/l Pb standard solution in 100 ml of distilled water, in a volumetric flask. Add also 20 μ l of Cd concentrated standard solution if cadmium has to be analysed in the same scanning.

Alternative supporting electrolytes

HCl or KCl or NaCl solution from 0.1 up to 1 M

0.1 M Acetate buffer pH 4.5 or 0.1 M citrate buffer a pH 3

0.1 M Tartrate buffer H 9 (when zinc has to be analysed in the same solution)

Avoid any addition of HNO_3 to the sample, because this acid could raise up the analytical peak. Use HCl instead of HNO_3 . If HNO_3 has to be used, dry the sample solution first and after add the supporting electrolyte.

Analytical report

Analysis: tap water

Sample Concentration = 4.40 $\mu\text{g/l}$

Method: 3 additions

Volumes Table

Solvent Volume	0 (ml)
Supporting Sol.	0.01 (ml)
Sample Volume	10 (ml)
Standard Conc.	1000 ($\mu\text{g/l}$)

Height Table

#	Peak Pot.	Height
0	-318.4	295.0 nA
1	-322.8	882.4 nA
2	-322.9	1.529 μA
3	-323.8	2.130 μA

Regression Data

#	Add.Conc.	Height x dilution	
0	0 $\mu\text{g/l}$	295.3 nA	$y = ax + b$
1	10.0 "	892.1 nA	$a = 63.75 \text{ nA}^*/\mu\text{g}$
2	20.0 "	1.562 μA	$b = 280.4 \text{ nA}$
3	30.0 "	2.197 μA	$r^2 = .9995$

