

## Palladium

**Method: 0.1 M Dimethylglyoxime in acetate buffer, pH 5.15**

**Function: Differential Pulse Stripping Voltammetry (DPS/a)**

Start Potential (mV)	-200
End Potential (mV)	-1000
Current range	1,024 $\mu$ A
Scan Speed (mV/s)	50
Deposition time (s)	120
Deposition Pot. (mV)	-200
Number of cycles	2
Delay before sweep (s)	5
Purge and stir time (s)	300
Stirring speed (rpm)	300
Drop Size (a.u.)	60

### Palladium concentrated standard Solution (1 g/l)

Dissolve 0.1 g of pure Pd in 5 ml of aqua regia (37% HCl + 65% HNO<sub>3</sub>, 3+1, v/v). Bring to dryness. Add of 37% HCl and 25 ml of distilled water. Heat gently, until salts are completely dissolved. Dilute to 100 ml with distilled water, in a volumetric flask.

### Supporting electrolyte

#### 1- 1 M acetate buffer, pH 5.15

Dissolve 13.6 g of CH<sub>3</sub>COONa · 3H<sub>2</sub>O (or 8.2 g of anhydrous CH<sub>3</sub>COONa) in 50 ml of distilled water, add 20% NaOH till pH 5.15. Dilute to 100 ml with distilled water, in a volumetric flask.

#### 2- 0.1 % Dimethylglyoxime in ethanol (p/v)

Dissolve 0.1 g of Dimethylglyoxime in 100 ml of ethanol. Prepare a fresh solution at the moment of the analysis.

### Procedure

Add to 10 ml of sample, 1 ml of buffer and 0.1 ml of dimethylglyoxime solution. Check that pH is 5.15. If necessary, correct by using NaOH or CH<sub>3</sub>COOH.

### Diluted standard solution (1 mg/l)

Dilute 1 + 999 the concentrated standard solution of Pd in distilled water. Prepare the solution at the moment of the analysis

### Working standard solution (100 $\mu$ g/l)

In a 50 ml volumetric flask, add 5 ml of buffer, 0.5 ml of dimethylglyoxime solution and 5 ml of diluted standard solution. Bring to volume with distilled water. Wait for 30 min before the use. Prepare the solution at the moment of the analysis.

### **Palladium in airborne**

#### **Procedure**

Sample the powder in the air using a cellulose filter, as described in the specific procedure for the determination of powder in air. Fold the filter and place it into the polarographic cell.

Add 2 ml of 65% HNO<sub>3</sub> and 2 ml of 40% H<sub>2</sub>O<sub>2</sub>. Let stand overnight.

Bring to dryness on a sand bath.

Add 1 ml of 65% HNO<sub>3</sub> and 2 ml of 40% H<sub>2</sub>O<sub>2</sub> and bring to dryness again.

Repeat the treatment until residue is white (not black, nor brown, nor yellow!)

Add 10 ml of distilled water to residue, 1 ml of buffer and 0.1 ml dimethylglyoxime solution. Check that pH is 5.15. If necessary, add NaOH or CH<sub>3</sub>COOH.

Alternatively, use a microwave digester but bring to dryness the residue.

#### **Warning**

Avoid using PTFE filters because the solution, after the boiling with concentrated HCl, cannot easily be digested.

## Analytical Report

Analysis: Filter n. 3 (143 mg of sampled powder  
in 288 l of air)

Sample Concentration = 0.64 mg Pd/kg of powder  
= 92 ng Pd on the filter  
= 0.32  $\mu\text{g}/\text{m}^3$  in the sampled air

Method: 4 additions

### Volumes Table

Solvent Volume	0 (ml)
Supporting Sol.	11.1 (ml)
Sample Weight	0.143 (g)
Standard Conc.	0.1 (mg/l)

### Height Table

#	Peak Pot.	Height
0	-593	145.0 nA
1	-610.7	428.5 nA
2	-606.8	622.0 nA
3	-616.2	749.4 nA
4	-616.2	881.9 nA

### Regression Data

#	Add. Conc.	Height x dilution
0	0 mg/kg	11.26 $\mu\text{A}$
1	1.40 "	39.26 $\mu\text{A}$
2	2.80 "	65.68 $\mu\text{A}$
3	4.20 "	89.62 $\mu\text{A}$
4	5.59 "	117.8 $\mu\text{A}$

$$y = ax + b$$

$$a = 18.84 \mu\text{A} \cdot \text{kg}/\text{mg}$$

$$b = 12.03 \mu\text{A}$$

$$r^2 = .9993$$

