

Sulphides

Function: Differential Pulse Voltammetry (DPV/a)

| | |
|-------------------------|-------|
| Start Potential (mV) | -200 |
| End Potential (mV) | -1000 |
| Current range | 4,096 |
| Scan Speed (mV/s) | 20 |
| Number of cycles | 3 |
| Delay before sweep (s) | 5 |
| Purge and stir time (s) | 300 |
| Stirring speed (rpm) | 300 |
| Drop Size (a.u.) | 60 |

Sulphide concentrated standard solution (1 g/l)

Dissolve 7.4901 g of $\text{Na}_2\text{S} \cdot 9\text{H}_2\text{O}$ in 1 l of 0.1 M NaOH, in a volumetric flask. Prepare the solution at the moment of the analysis. ($\text{MM}_{\text{Na}_2\text{S}} = 240.2$ $\text{MM}_{\text{S}} = 32$)

Supporting electrolyte

1 M NaOH solution

Dissolve 40 g of NaOH in 1 l of distilled water.

Procedure

Pour 10 ml of sample in the cell. Add 1 ml of supporting electrolyte.

Working standard solution (10 mg/l)

Dilute 1+99 the concentrated standard solution with distilled water. Prepare the solution at the moment of the analysis

Warnings

- Store samples in hermetically sealed bottles without bubble air inside. Analyse these samples as soon as possible
- For a better treatment, add 4 g of NaOH for 1 l to the sample, at the moment of sampling and analyse the solution as soon as possible.

Analytical report

Analysis: Anaerobic water plant
 Sample Concentration = 14.9 $\mu\text{g/l}$
 Method: 5 additions

Volumes Table

| | |
|-----------------|---------------------------|
| Solvent Volume | 0 (ml) |
| Supporting Sol. | 1 (ml) |
| Sample Volume | 10 (ml) |
| Standard Conc. | 10000 ($\mu\text{g/l}$) |

Height Table

| # | Peak Pot. | Height |
|---|-----------|---------------------|
| 0 | -637.5 | 1.316 μA |
| 1 | -657.5 | 8.259 μA |
| 2 | -663.8 | 15.50 μA |
| 3 | -671.8 | 22.55 μA |
| 4 | -676.3 | 29.37 μA |
| 5 | -681.2 | 36.57 μA |

Regression Data

| # | Add. Conc. | Height x dilution |
|---|-------------------|---------------------|
| 0 | 0 $\mu\text{g/l}$ | 1.448 μA |
| 1 | 100 " | 9.168 μA |
| 2 | 200 " | 17.37 μA |
| 3 | 300 " | 25.49 μA |
| 4 | 400 " | 33.49 μA |
| 5 | 500 " | 42.06 μA |

$$y = ax + b$$

$$a = 81.19 \text{ nA}^*/\mu\text{g}$$

$$b = 1.207 \mu\text{A}$$

$$r^2 = .9998$$

