

Sulphur dioxide - Sulphides

Function: Differential Pulse Voltammetry (DPV/a)

Start Potential (mV)	-300
End Potential (mV)	-900
Current range	1,024 μ A
Scan Speed (mV/s)	20
Number of cycles	3
Delay before sweep (s)	5
Purge and stir time (s)	300
Stirring speed (rpm)	300
Drop Size (a.u.)	60

Concentrated standard solution (1 g/l) of sulphur dioxide

Dissolve 1.9688 g of anhydrous Na_2SO_3 in 1 l of distilled and boiled water, in volumetric flask. Prepare the solution at the moment of the use. ($\text{MM}_{\text{Na}_2\text{SO}_3} = 126.04$ $\text{MM}_{\text{SO}_2} = 64$)

Supporting electrolyte

1 M Acetate buffer, pH 4.6

Dissolve 82 g of CH_3COONa (or 136 g of $\text{CH}_3\text{COONa} \cdot 3\text{H}_2\text{O}$) in 800 ml of distilled water. Add 57.5 ml of glacial CH_3COOH . Adjust pH and bring to the mark in a 1 l volumetric flask

Procedure

Deaerate 10 ml of supporting electrolyte and add from 1 to 5 ml of sample.

Working standard solution (100 mg/l)

Dilute 1+9 the concentrated standard solution with distilled and boiled water. Prepare the solution at the moment of the analysis

Warnings

- Store acidic samples in hermetically sealed bottles without bubble air inside. Analyse these samples as soon as possible
- If sulphides are present (>0.2 mg/l), start the scanning from -400 , -450 mV.
- Dilute concentrated samples using boiled distilled water.

Analytical report

Analysis: wine Barbera

Sample Concentration = 164 mg/l

Method: 5 additions

Volumes Table

Solvent Volume	2.6 (ml)
Supporting Sol.	10 (ml)
Sample Volume	0.5 (ml)
Standard Conc.	100 (mg/l)

Heights Table

#	Peak Pot.	Height
0	-597	1.100 μA
1	-594	1.612 μA
2	-592.5	2.185 μA
3	-592.5	2.684 μA
4	-589.5	3.119 μA
5	-588.6	3.643 μA

Regression Data

#	Add. Conc.	Height x dilution
0	0 mg/l	28.83 μA
1	100 "	43.86 μA
2	200 "	61.62 μA
3	300 "	78.38 μA
4	400 "	94.22 μA
5	500 "	113.7 μA

