


Karl Fischer Titration Electrodes

Table of Contents

1	880/TPG/12J.....	1
2	880/DPE/12J.....	1

1 880/TPG/12J

Double platinum electrode with independent metallic endings, comes with S7 connector only. Conical tapered joint and 12mm body diameter to fit AMEL's Karl Fischer titrators.



2 880/DPE/12J

Double platinum electrode with independent metallic endings, comes with S7 connector only. Conical tapered joint and 9mm final body diameter to fit common Karl Fischer titrators.

